

We simplify complex processes in



## Digital Wellness Assessment

Wellness Score based on

Lifestyle Markers

Autonomic Nervous System Aging

Vascular Aging

Bodyk.net

[bodyknows@bodyk.net](mailto:bodyknows@bodyk.net)

Miami, FL - USA





## BodyK DIGITAL ASSESSMENT

---

### WHAT IS **BodyK** ?

**BodyK** is a revolutionary cutting-edge technology for functional body assessments with the fastest, non-invasive, and most comprehensive, validated wellness screening, positioning your practice at the forefront of technology. We get to the root-cause of your patients' symptoms in just 7 minutes.


**BodyK IS THE WAY TO HELP YOUR PATIENTS REACH THEIR OPTIMAL HEALTH,  
WHILE INCREASING YOUR PRACTICE'S MONTHLY REVENUE.**



# INTEGRATED TECHNOLOGIES



GALVANIC SKIN RESPONSE  
BIOIMPEDANCE ANALYSIS




OXYMETER



PHOTOPLETHYSMOGRAPHY



BLUETOOTH BLOOD PRESSURE



FDA CLEARED AS A GENERAL WELLNESS DEVICE

PATENTED TECHNOLOGIES

Manufacturer Certified






# WMS TESTS METHOD AND PROCEDURE

## PROCEDURE:

1. Patient is sitting and relax (at least for 5 min)
2. Set up patient – update software with patient info.
3. Place the disposable electrodes on the feet and forehead, then place cuff on the left arm, and finally place pulse oximeter on the finger of the right arm.
4. Start the exam and follow the instructions on the software.


## 1- Setup patient



## 2-RECORDING



## 3-RESULTS





## BODYK SYSTEM

---

OUTPUT DATA



# BODYK OUTPUT DATA

**Visit Date:** 11/11/2023      **Visit Time:** 14:58  
**Client Name:** [REDACTED]  
**Gender:** Male      **Height:** 5' 11"      **BMI:** 33.2  
**DOB:** 8/27/1986      **Weight:** 237 Lbs  
**Age:** 37      **Vascular Age:** 42

### WELLNESS SCORE

**62%**

SCORE 0% - 100% \*

100%	EXCELLENT
90%	GOOD
80%	BORDERLINE
60%	POOR
<= 40%	VERY POOR

## Client Wellness Profile

### WELLNESS ASSESSMENTS

 <b>INSULIN RESISTANCE</b> <b>DETECTED</b>	 <b>Activity Level</b> <b>Moderate</b> <b>Heart Rate</b> <b>86 bpm</b> <b>BMI</b> <b>33.2</b> <b>SpO2%</b> <b>93</b>	 <b>FAT MASS PERCENT</b> <b>ABNORMAL</b>
 <b>INFLAMMATION</b> <b>BORDERLINE</b>		 <b>BLOOD PRESSURE</b> <b>PRE-HYPERTENSION</b>
 <b>IMMUNE RESPONSE</b> <b>NORMAL</b>		 <b>OXIDATIVE STRESS</b> <b>BORDERLINE</b>
 <b>FITNESS SCORE</b> <b>GOOD</b>		 <b>MICROCIRCULATION</b> <b>BORDERLINE</b>
 <b>STRESS LEVEL</b> <b>BORDERLINE</b>		 <b>DIGESTIVE INFLAMMATION</b> <b>NORMAL</b>

**ASSESSMENTS COLOR CODE**

█ NORMAL RANGE   
 █ INITIAL STEP 1 / ACCEPTABLE   
 █ INITIAL STEP 2 / BORDERLINE   
 █ DEFINED   
 █ SEVERE



# BODYK OUTPUT DATA

The BodyK outputs are an overview of markers related to a Wellness risk calculated as a score.

## Results

### WELLNESS SCORING

### WELLNESS ASSESSMENTS

### ASSESSMENT MEANING



17 % LIFESTYLE SCORE \*\*

22 % AUTONOMIC SCORE \*\*\*

23 % VASCULAR AGING SCORE \*\*\*

### Main Markers Overview

#### LIFESTYLE

FITNESS BIOMARKERS	VISIT 1
CARDIOVASCULAR FITNESS <i>Total Power &gt;= 780 (ms2)</i>	809
PERFORMANCE CAPACITY <i>SDANN Age dependent (ms)</i>	31
EXERCISE RECOVERY INDEX <i>rMSSD Age dependant (ms)</i>	26
BLOOD OXYGEN SATURATION <i>Spo2 &gt;= 94 (%)</i>	93
MUSCLE AND BONE % <i>Age and gender dependent</i>	27
<b>DIET</b>	
FAT MASS % <i>Age and gender dependent</i>	41
BODY MASS INDEX <i>Age / Gender dependent.</i>	33
INTERSTITIAL PH REGULATION <i>Delta C &gt;= 40 and &lt;= 79 (mSi)</i>	49
<b>EMOTIONAL</b>	
MENTAL STRESS <i>LF / HF &lt; 2 (ratio)</i>	2.3
RECOVERY CAPACITY <i>HF &gt;= 22 (%)</i>	19
SEROTONIN <i>Head Nitric Oxide &gt; 256 (mV)</i>	1178
DOPAMINE <i>Head Sweat Peak &lt;= 768 (mV)</i>	1088

#### AUTONOMIC NERVOUS SYSTEM

SYMPATHETIC RESPONSE TESTS	VISIT 1
PRESSURE/ HEART RATE REGULATION <i>Delta LF &gt;= 1.2 (ratio)</i>	1.5
NORADRENERGIC RESPONSE <i>SPRS &lt; 10 (mmHg)</i>	36
<b>PARASYMPATHETIC RESPONSE TEST</b>	
HEART RATE RESPONSE TO CHALLENGE <i>E/I R. Age/Gender dependent</i>	1.25
HEART RATE RESPONSE TO STANDING <i>K3015 R. Age/Gender dependent</i>	1.24
<b>GSR TESTS</b>	
MICROCIRCULATION <i>Foot Nitric Oxide &gt;= 768 (mV)</i>	614
OXIDATIVE STRESS <i>Foot Sweat Peak &gt;= 832 (mV)</i>	627
THYROID FUNCTION <i>Upper voltage &gt;= 512 (mV)</i>	1133
DIGESTIVE INFLAMMATION <i>Middle voltage &gt;= 512 and &lt;= 768 (mV)</i>	614

#### VASCULAR FUNCTION

BLOOD PRESSURE	VISIT 1
SYSTOLIC PRESSURE <i>&lt; 130 (mmHg)</i>	127
DIASTOLIC PRESSURE <i>&lt; 85 (mmHg)</i>	85
STROKE VOLUME <i>&gt;= 60 and &lt;= 100 (ml/beat)</i>	68.8
CARDIAC OUTPUT <i>&gt;= 4 and &lt;= 8 (ml/min)</i>	5.9
<b>ENDOTHELIAL FUNCTION</b>	
INFLAMMATION <i>Stress Index &lt; 180 (%)</i>	228
IMMUNE RESPONSE <i>PTG Index &gt;= 40 (Vs)</i>	62.6
LIPID PROFILE <i>PTGAI &lt;= 0.45 (ratio)</i>	0.17
VASCULAR TONE <i>SDPTG &lt;= 0.42 (ratio)</i>	0.21
INSULIN RESISTANCE <i>PTG TP &lt;= 406 (ms2)</i>	540
COAGULATION <i>PTGVLF Index &lt;= 32 (Vs / mS)</i>	85

(*)	(**)	(***)	
100%	40%	30%	EXCELLENT
90%	≥ 37%	≥ 28%	GOOD
80%	≥ 32%	≥ 25%	BORDERLINE
60%	≥ 24%	≥ 18%	POOR
≤ 40%	≤ 16%	≤ 12%	VERY POOR

#### COLOR CODE

■ NORMAL RANGE   
 ■ BORDERLINE   
 ■ ABNORMAL

CURRENT VISIT 1

11/11/2023 14:58

MARKERS OVERVIEW



# BODYK OUTPUT DATA

## Marker overview and comparison with another exam

## Score trends

### Main Markers Overview

#### LIFESTYLE

FITNESS BIOMARKERS	VISIT 1	VISIT 2
CARDIOVASCULAR FITNESS <i>Total Power &gt;= 780 (ms2)</i>	1835	1264
PERFORMANCE CAPACITY <i>SDANN Age dependent. (ms)</i>	53	42
EXERCISE RECOVERY INDEX <i>rMSSD Age dependant (ms)</i>	64	43
BLOOD OXYGEN SATURATION <i>SpO2 &gt;= 94 (%)</i>	96	97
MUSCLE AND BONE % <i>Age and gender dependent</i>	36	36
<b>DIET</b>		
FAT MASS % <i>Age and gender dependent</i>	18	19
BODY MASS INDEX <i>Age / Gender dependent</i>	22.5	23.1
INTERSTITIAL PH REGULATION <i>Delta C &gt;= 40 and &lt;= 79 (mSI)</i>	94	91
<b>EMOTIONAL</b>		
MENTAL STRESS <i>LF / HF &lt; 2 (ratio)</i>	0.9	3.4
RECOVERY CAPACITY <i>HF &gt;= 22 (%)</i>	36	20
SEROTONIN <i>Head Nitric Oxide &gt; 256 (mV)</i>	1242	1242
DOPAMINE <i>Head Sweat Peak &lt;= 768 (mV)</i>	909	1114

#### AUTONOMIC NERVOUS SYSTEM

SYMPATHETIC RESPONSE TESTS	VISIT 1	VISIT 2
PRESSURE/ HEART RATE REGULATION <i>Delta LF &gt;= 1.2 (ratio)</i>	1.5	0.8
NORADRENERGIC RESPONSE <i>SPRS &lt; 10 (mmHg)</i>	8	16
<b>PARASYMPATHETIC RESPONSE TEST</b>		
HEART RATE RESPONSE TO CHALLENGE <i>E/I R. Age/Gender dependent</i>	1.16	1.26
HEART RATE RESPONSE TO STANDING <i>K3015 R. Age/Gender dependent</i>	1.30	1.24
<b>GSR TESTS</b>		
MICROCIRCULATION <i>Foot Nitric Oxide &gt;= 768 (mV)</i>	1114	1037
OXIDATIVE STRESS <i>Foot Sweat Peak &gt;= 632 (mV)</i>	1280	1280
THYROID FUNCTION <i>Upper voltage &gt;= 512 (mV)</i>	1075	1178
DIGESTIVE INFLAMMATION <i>Middle voltage &gt;= 512 and &lt;= 768 (mV)</i>	1190	1152

#### VASCULAR FUNCTION

BLOOD PRESSURE	VISIT 1	VISIT 2
SYSTOLIC PRESSURE <i>&lt; 130 (mmHg)</i>	134	132
DIASTOLIC PRESSURE <i>&lt; 85 (mmHg)</i>	74	75
STROKE VOLUME <i>&gt;= 60 and &lt;= 100 (ml/beat)</i>	120	114
CARDIAC OUTPUT <i>&gt;= 4 and &lt;= 8 (ml/min)</i>	6.4	8
<b>ENDOTHELIAL FUNCTION</b>		
INFLAMMATION <i>Stress Index &lt; 180 (%)</i>	54	95
IMMUNE RESPONSE <i>PTG Index &gt;= 40 (Vs)</i>	44.8	43.4
LIPID PROFILE <i>PI(GAI) &lt;= 0.45 (ratio)</i>	0.34	0.35
VASCULAR TONE <i>SDPTG &lt;= 0.42 (ratio)</i>	0.78	0.37
INSULIN RESISTANCE <i>PTG TP &lt;= 406 (ms2)</i>	269	372
COAGULATION <i>PTGVLf Index &lt;= 32 (Vs / mS)</i>	12	32


COLOR CODE  
■ NORMAL RANGE ■ BORDERLINE ■ ABNORMAL

CURRENT VISIT 1 9/25/2021 12:10  
 COMPARE VISIT 2 5/6/2021 16:41



## Charts


### BODY COMPOSITION CHART




### FITNESS CHART




### MICROCIRCULATION AND OXIDATIVE STRESS RISK CHART




### ENDOTHELIAL DYSFUNCTION RISK CHART



### CARDIAC AUTONOMIC NEUROPATHY (CAN) RISK CHART



### HYPERTENSION RISK CHART










# BODYK SYSTEM

---

GALVANIC SKIN RESPONSE  
TECHNOLOGY






# SWEATC : PROCESS OF MEASUREMENT



**NO Sweat Peak Marker of microvascular disorders and PH**

Current generate by the hardware




**Sweat Peak marker of Oxidative stress damaging the sympathetic cholinergic fibers**



## SWEATC CLINICAL EVALUATION

ROC Curve FOR NO Sweat Peak



### NO SWEAT PEAK MARKER OF MICROVASCULAR DISORDERS :

**NO Sweat Peak inversely correlated with :**

**Lab tests:**

**BUN** (Spearman  $\rho=-0.41$ ,  $p<0.0001$ )

**Homocysteine** (Spearman  $\rho=-0.44$ ,  $p<0.0001$ )

**Fibrinogen** (Spearman  $\rho=-0.41$ ,  $p<0.0001$ )

**Neuropathy (vascular disease)**

**CAN score** (Spearman  $\rho=-0.68$ ,  $p<0.0001$ )

**HRV Total Power** (Spearman  $\rho=-0.57$ ,  $p<0.0001$ )

**Marker of endothelial function:**

**PTGi** (Spearman  $\rho=0.53$   $p<0.0001$ )

**NO Sweat Peak had a sensitivity of 88% and specificity of 68% to detect retinopathy ( microvascular disease)**

### ISWEAT PEAK MARKER OF PERIPHERAL NEUROPATHY

**iSweat Peak inversely correlated with the severity of symptoms on the peripheral neuropathy scale (Spearman  $\rho=-0.56$ ,  $p<0.0001$ ).**

University of Miami study . J. Lewis et al. 2017



# BODYK SYSTEM


---

PHOTOPLETHYSMOGRAPHY  
TECHNOLOGY



# PHOTOPLETHYSMOGRAPHY

Incident Light (Red/IR)




PTG :Arterial vessel and changes in blood volume with cardiac cycle

Venous vessel not affected




# PHOTOPLETHYSMOGRAPHY (PTG) RECORDS AND ANALYSIS


First derivative PTG Analysis



Original PTG Analysis





Second derivative PTG Analysis






# Peripheral vascular tone related to age




## Pulse Waveform of a 47 year old person

With age, the arteries stiffen. Pulse wave velocity increases and the reflected wave now travels faster and gives rise to a "shouldering" at the primary wave, and occurs during systole. This effect prolongs the systolic cycle increasing the workload and oxygen requirement of the heart muscles. The poor notch also compromises the flow of the coronaries.






# PHOTOPLETHYSMOGRAPHY ANALYSIS BACKGROUND




A  
N  
A  
L  
Y  
S  
I  
S

## 1 ORIGINAL PHOTOPLETHYSMOGRAPHY (PTG)




1

## 2 FIRST DERIVATIVE AND RR INTERVALS: TIME IN MS BETWEEN EACH HEARTBEAT




2

## 3 PTG SECOND DERIVATIVE



3

## 4 PATENTED FIRST DERIVATIVE SPECTRAL ANALYSIS




4





# RR INTERVALS ANALYSIS AND BEAT TO BEAT BLOOD PRESSURE ANALYSIS

## HEART RATE VARIABILITY ANALYSIS AT REST



Self body Regulation Assessment related to the fitness level


## HEART RATE VARIABILITY AND BLOOD PRESSURE ANALYSIS DURING CHALLENGES



Diagnostic of Cardiac Autonomic Neuropathy (CAN)



# PATENTED PTG SPECTRAL ANALYSIS AND MARKERS



**PTG INDEX (PTGI) =  
MARKER OF ENDOTHELIAL FUNCTION**  
It is the sum of amplitude between the  
Peaks of the Spectral analysis.

**STRESS INDEX =  
MARKER OF INFLAMMATION**  
It is the amplitude of the VLF  
peak of the Spectral analysis.

**PTG TOTAL POWER (PTG-TP) =  
MARKER OF INSULIN RESISTANCE**  
It is the area covering the 3 frequencies  
of the PTG spectral Analysis



## BODYK SYSTEM


---

PATENTED PTG SPECTRAL ANALYSIS AND  
CLINICAL OUTCOMES



# Diabetes Pathogenesis


## DIABETES PROCESS STARTED 10 YEARS BEFORE DIAGNOSIS





## Insulin resistance and Health risk

1. Insulin resistance is the higher risk to be obese
2. 83% of patients with Insulin resistance will go in Diabetes
3. Insulin resistance is higher risk factor for cardiovascular disease : greater than smoking, cholesterol and aging



Bonora E, et al. *Diabetes Care* 2002; 25:1135–1141

Insulin resistance could be reversed by 58% following lifestyle changes



## CLINICAL STUDIES – INSULIN RESISTANCE

Aglecio Luiz De Souza, Gisele Almeida Batista, Sarah Monte Alegre  
*Journal of Diabetes and its complications.*



### ASSESSMENT OF INSULIN SENSITIVITY BY THE HYPERINSULINEMIC EUGLYCEMIC CLAMP: COMPARISON WITH THE SPECTRAL ANALYSIS OF PHOTOPLETHYSMOGRAPHY

#### RESULTS

Correlation between insulin sensitivity (M-value) and PTG-TP ( $r = -0.64$ ,  $p < 0.0001$ ). PTG-TP had a sensitivity = 95.7%, specificity = 84.4% and the area under the ROC curve (AUC) = 0.929 for identifying insulin resistance. ( $p < 0.0001$ ).

#### CONCLUSION


The use of the PTG-TP marker measured from the PTG spectral analysis is a useful tool in screening and follow up of Insulin Resistance (IR)





## INSULIN RESISTANCE DETECTION PERFORMANCE

### Insulin resistance detection marker (PTG-TP) vs. Gold standard and other alternative blood tests and algorithms



*DE SOUZA ET AL. Assessment of insulin sensitivity by the hyperinsulinemic euglycemic clamp: Comparison with the spectral analysis of Photoplethysmography. Journal of Diabetes and its complications. Volume 31, (2017) 128–133*



## CLINICAL STUDIES – ENDOTHELIAL FUNCTION

Gandhi PG, Rao GHR.

### THE SPECTRAL ANALYSIS OF THE PHOTOPLETHYSMOGRAPHY TO EVALUATE AN INDEPENDENT CARDIOVASCULAR RISK. INT J GEN MED.

#### RESULTS



PTGi had a sensitivity of 86.5% and specificity of 87.8%, ( $P=0.0001$ ; area under the receiver operating characteristic curve =0.967) to detect endothelial dysfunction in CAD population.

Patients that underwent CABG ( $n=18$ ) had a higher PTGi value compared with the CAD without CABG surgery patients ( $n=47$ ).

#### CONCLUSION

The spectral analysis techniques used on the photoplethysmogram, as outlined in this study, could be useful when used alongside conventional known cardiovascular disease risk markers.

#### RECEIVER OPERATING CHARACTERISTIC CURVES








*UM, Lewis et al. 2017*

## A NEW METHOD OF SUDOMOTOR FUNCTION MEASUREMENT TO DETECT MICROVASCULAR DISEASE AND SWEAT GLAND NERVE OR UNMYELINATED C FIBER DYSFUNCTION IN ADULTS WITH RETINOPATHY.

### RESULTS

The marker Stress Index correlated with C-Reactive Protein (Spearman  $\rho=0.40$ ,  $p<0.0001$ )





# BODYK SYSTEM

---

## MODELING REFERENCES



## PH , OXIDATIVE STRESS AND EMOTIONAL MARKERS

**INTERSTITIAL PH** is estimated from the delta of conductance of the left and right plates of the SweatC device. The scientific support is described in the following publication: Maarek A. Electro interstitial scan system: assessment of 10 years of research and development. Medical Devices: Evidence and Research. March 2012 Volume 2012:5 Pages 23 - 30.

---

**OXIDATIVE STRESS** is estimated from the conductance average of the Sweat Peak of the SweatC device. The scientific support is described in the following publication: Oxidative Stress in the Vincent A.M, Russell J.W, Low P, Feldman E. L. Pathogenesis of Diabetic Neuropathy Endocrine Reviews, Volume 25, Issue 4, 1 August 2004, Pages 612–628,

---

**BRAIN SEROTONIN** is estimated from the conductance average of the left and right forehead of the SweatC device. The scientific support is described in the following publication: VG ALEXEEV VG, LV KUZNECOV. Bioimpedance in monitoring of effects of selective serotonin reuptake inhibitor treatment. Psychology Research and Behavior Management. p.81-86. Volume 4.2011.

---

**MENTAL STRESS** is estimated from the marker LF/HF of the HRV analysis. The scientific support is described in the following publication: Taelman J, Vandepuut S , Spaepen A, Van Huffel S.4th European Conference of the International Federation for Medical and Biological Engineering pp 1366-1369 | Cite asInfluence of Mental Stress on Heart Rate and Heart Rate Variability.

---

**BRAIN DOPAMINE** is estimated from the conductance average of the left and right forehead of the SweatC device. The scientific support is described in the following publication: F. CAUDAL. New marker using bioimpedance technology in screening for attention deficit/hyperactivity disorder (ADHD) in children as an adjunct to conventional diagnostic methods. Psychology Research and Behavior Management » Volume 4, P.113-117.2011



## HYPOTHYROIDISM AND DIGESTIVE SYSTEM

**HYPOTHYROIDISM** is evaluated from the autonomic nervous system markers and adjusted with some side effects of the hypofunction such as , blood pressure and arterial stiffness . The scientific support is described in the following publication: Mahajan AS, Lal R, Dhanwal DK, Jain AK, and Chowdhury V. Evaluation of autonomic functions in subclinical hypothyroid and hypothyroid patients. Indian J Endocrinol Metab. 2013 May-Jun; 17(3): 460–464.

---

**DIGESTIVE SYSTEM** is evaluated from the autonomic nervous system markers. The scientific support is described in the following publication: Browning KN, Travagli RA. Central nervous system control of gastrointestinal motility and secretion and modulation of gastrointestinal functions. Compr Physiol. 2014 Oct;4(4):1339-68.

---

**HEPATIC GLYCOLYSIS** is evaluated from the autonomic nervous system markers. The scientific support is described in the following publication: Nonogaki K. New insights into sympathetic regulation of glucose and fat metabolism. Diabetologia. 2000 May;43(5):533-49



# BODYK SYSTEM

---

TECHNOLOGIES SUPPORTIVE STUDIES



# MAIN CLINICAL STUDIES

- [A.MAAREK.](#)  
Endothelial and Autonomic Assessments in the Management of COVID-19.  
Open Journal of Preventive Medicine Vol.11 No.12, December 23, 2021
- [LOPEZ ET AL.](#)  
Relationship Between Vitamin D Status, Autonomic nervous system activity and Cardiometabolic Risk.  
Open Journal of internal medicine, 2017,7,37-51
- [LEWIS ET AL.](#)  
New method of sudomotor function measurement to detect microvascular disease and sweat gland nerve or unmyelinated C fiber dysfunction in adults with retinopathy. Journal of Diabetes & Metabolic Disorders (2017) 16:26.
- [DE SOUZA ET AL.](#)  
Assessment of insulin sensitivity by the hyperinsulinemic euglycemic clamp: Comparison with the spectral analysis of photoplethysmography. Journal of Diabetes and its complications. Volume 31, (2017) 128–133.
- [LEWIS ET AL.](#)  
A cross-sectional assessment to detect Type 2 Diabetes with endothelial and autonomic nervous system markers using a novel system.  
Journal of Diabetes & metabolic disorders, 2014,13:118
- [GOBATO RC, SEIXAS CHAVES DF, CHAIM EA.](#)  
Micronutrient and physiologic parameters before and 6 months after RYGB. Surg Obese Relat Dis. 2014 Sep-Oct;10(5):944-51.
- [GANDHI PG, RAO GHR.](#)  
The spectral analysis of the photoplethysmography to evaluate an independent cardiovascular risk.  
Int J Gen Med. 2014 Dec 9;7:539-47.
- [GANDHI PG, RAO GHR.](#)  
Detection of neuropathy using a sudomotor test in type 2 Diabetes. Degenerative Neurological and Neuromuscular Disease. January 2015, 5:1-7.
- [MAAREK A A., GANDHI PG AND RAO GHR.](#)  
Identifying Autonomic Neuropathy and Endothelial Dysfunction in Type II Diabetic Patients. EC Neurology 2.2(2015): 63-78.
- [GUNDU HR RAO, PRATIKSHA G GANDHI AND VINEETA SHARMA.](#)  
Clinical Complications of Type-2 Diabetes Mellitus in South Asian and Chinese Populations: An Overview. J Diabetes Metab. Volume 5 • Issue 8.
- [GUNDU H R RAO, PRATIKSHA G GANDHI.](#)  
Need for a Noninvasive Diagnostic Platform for Early Detection and Management of Cardiometabolic Disorders. Journal of Clinical and Preventive Cardiology, July 2014 | Number 3.
- [GUNDU HR RAO.](#)  
Flow Velocity, Fluid Dynamics and Vascular Pathophysiology. The scientific pages of Heart.2016,1.001,Volume 1.Issue1.
- [MATJAZ Š, GERSAK G, MILLASSEAU SC, MEZA M, KOSIR A.](#)  
Detection of peripheral arterial disease with an improved automated device: comparison of a new oscillometric device and the standard Doppler method. Vascular Health and Risk Management.2016. Volume 12.
- [RAO GHR.](#)  
Non-Traditional Approaches to Diagnosis and Management of Type-2 Diabetes Mellitus: Point of View. J Diabetes Metab (2015) 6: 489.



## R & D PUBLICATIONS

- [MAAREK A.](#)

Electro interstitial scan system: assessment of 10 years of research and development.

Medical Devices: Evidence and Research.

[March 2012 Volume 2012:5, Pages 23 – 30.](#)

- [F. CAUDAL](#)

New marker using bioimpedance technology in screening for attention deficit/hyperactivity disorder (ADHD) in children as an adjunct to conventional diagnostic methods. Psychology Research and Behavior Management » Volume 4, P.113-117.2011

- [VG ALEXEEV VG, LV KUZNECOV.](#)

Bioimpedance in monitoring of effects of selective serotonin reuptake inhibitor treatment. Psychology Research and Behavior Management. p.81-86. Volume 4.2011.

## PEER REVIEWS

- [Photoplethysmography and its application in clinical physiological measurement](#)
- [Diabetes Autonomic Neuropathy](#)
- [Diabetic Neuropathy: Diagnostic Methods](#)
- [Determination of age-related increases in large artery stiffness by digital pulse contour analysis](#)
- [Role of Endothelial Dysfunction in Atherosclerosis](#)
- [Impact of Reduced Heart Rate Variability on Risk for Cardiac Events](#)
- [Diabetes, Glucose, Insulin, and Heart Rate Variability](#)
- [Diabetic Neuropathies: The Nerve Damage of Diabetes](#)
- [Diabetes, Glucose, Insulin, and Heart Rate Variability](#)
- [Independent Determinants of Second Derivative of the Finger Photoplethysmogram among Various Cardiovascular Risk Factors in Middle-Aged Men](#)
- [Assessment of Vascular Aging and Atherosclerosis in Hypertensive Subjects](#)
- [Diabetes: Sweat Response and Heart Rate Variability During Electrical Stimulation in Controls and People With Diabetes](#)
- [Diagnosis and Treatment of Pain in Small Fiber Neuropathy](#)
- [Fitness markers – Department of Analytics and Physiology](#)
- [HRV fitness markers](#)
- [Physical capacity with heart rate variability](#)